

# The Sercon CryoGas: for automated, high precision $\delta^2\text{H}$ and $\delta^{13}\text{C}$ analysis of gaseous and dissolved methane



Andrew C. Smith<sup>1\*</sup>, Steve Welsh<sup>2</sup>, Helen Atkinson<sup>2</sup> and Melanie Leng<sup>1</sup>

<sup>1</sup> National Environmental Isotope Facility, British Geological Survey, Nottingham, NG12 5GG

<sup>2</sup> Sercon Group, 3B, Crewe Trade Park, Gateway, Crewe CW1 6JT

\* email: andrews@bgs.ac.uk




## 1. Introduction:

- Methane is the **third most important greenhouse gas** after water vapour and  $\text{CO}_2$ , contributing 20% of all radiative forcing from the long lived green house gases.
- Isotopic analysis** of carbon ( $^{13}\text{C}/^{12}\text{C}$ ) and hydrogen ( $^2\text{H}/^1\text{H}$ ) in methane, can tell us by which process the methane has been formed. Enabling us to understand where methane pollution is originating from and mitigate against it.
- The new **Sercon CryoGas** (Figure 1) is an automated system for the high precision, rapid throughput analysis of both carbon and hydrogen isotopes in methane.




**Figure 1:** The Sercon CryoGas and HS2022 mass spectrometer, this is the first Sercon instrument built specifically for the analysis of C and H in methane.



## 2. Operation:


- Samples are introduced as a mixed gas, either atmosphere or an equilibrated headspace from water samples.
- In both hydrogen and carbon modes the CryoGas utilises **cryogenic trapping, chemical gas clean up and gas chromatography** to separate out air components ( $\text{CO}_2$ , N,  $\text{O}_2$ ) from methane (Figures 2 and 3).
- Pyrolysis is used for  $\delta^2\text{H}$  - $\text{CH}_4$  and combustion for  $\delta^{13}\text{C}$ - $\text{CH}_4$  analysis.



**Figure 2:** Schematics of the CryoGas system in both carbon (top) and hydrogen (bottom) modes. Showing a series of chemical and cryo traps as well as the configuration of GC columns/ vents for the separation of  $\text{CH}_4$  from air components.

## 3 Results:

- Instrument precision** on 150ml of air is  $<3\text{‰}$  for  $\delta^2\text{H}$  - $\text{CH}_4$  (Figure 4b) and  $<0.2\text{‰}$  for  $\delta^{13}\text{C}$ - $\text{CH}_4$ .
- The CryoGas can do **unattended runs** of either 21 samples of 150ml or  $>200$  12ml samples at higher concentrations.



**Figure 3a)** A typical  $\delta^2\text{H}$  - $\text{CH}_4$  in air output, with initial ref gas pulse followed by a  $>1$  nano amp sample peak height.  
**3b)** example data showing the CryoGas precision for 150ml of air and standards at 2000ppb  $\text{CH}_4$ .

